
AESTHETIC-PRESS
EFFICIENCY THROUGH SIMPLICITY

A Perfect Workflow by Technicians for Technicians

Product Catalog / Produktkatalog

AESTHETIC PRESS SEMINARS / WORKSHOPS

LECTURER / INSTRUCTOR BIOGRAPHY

Jörg Müller
MDT, President of Aesthetic-Press

Before receiving his degree as a master technician in 1993, Jörg was trained in Jan Langner's Dental Laboratory in Schwäbisch-Gmünd, where he specialized in tooth morphology function and aesthetics.

In 1993, working closely with Claude Sieber for several years on anterior aesthetics and photography, Jörg opened his first dental laboratory in Düsseldorf. Since then, he has been teaching courses on posterior wax-up, the function of the TMJ, and ceramics for dental technicians, worldwide.

In 1999 Jörg moved to San Francisco, California and established his dental laboratory in downtown San Francisco. While working closely with his dentists and patients, he developed the Aesthetic Press System. In 2005, Jörg founded Aesthetic-Press LLC with which he developed a system to standardize a workflow for consistent and reproducible high-end restorations for technicians to use on a day-to-day basis.

Today, Aesthetic Press is based in Düsseldorf, Germany and Palm Harbor, FL.

Seminars and workshops are available for continuing education credits.

Available courses:

- Introduction to AP and Pressing (5 CE credit hours)
- Aesthetic-Press Workshop: Mastering the Press Technique (15 CE credit hours)

Contact us for an updated list of available seminars and workshops

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81-84

INVESTMENT / EINBETTEN

Easy Vest Investment
 Easy Vest Investment PREMIUM (NEW)
 Plunger Maker / *Pressstempelformer*
 Investment Ring / *Muffel*

85-86

AQUALINE

Aqualine Smart
 Aqualine Mini
 Aqualine Stain Tray

87

DEKEMA OVENS

AUSTROMAT 624
 AUSTROMAT 654 press-i-dent
 AUSTROMAT 664
 AUSTROMAT 674
 AUSTROMAT μ SIC
 AUSTROMAT 624 i
 AUSTROMAT 654 i
 AUSTROMAT 664 i
 AUSTROMAT 674 i

88-89

ROTARY SYSTEMS

Diamond Burs
 Carbide Burs
 Stone Grinders
 Panther
 Diamond Discs
 Polishers
 CAD/CAM
 Kits

Simplicity is the ultimate sophistication!

Leonardo DaVinci

From Technician to Technician!

Jörg Müller, founder and developer of Aesthetic-Press:

During my career as a dental technician, the ultimate goal was to create a consistent, high quality standard for each case every technician produced in the lab. With the AP system, we've achieved a consistent, predictable and reproducible level of the highest quality. The simplicity of the integrated product line is one of its strongest points.

Next to our unique inventions, we've improved and streamlined the workflow into a smooth product line without clutter or more than a technician needs to be successful. In essence, we have made dental technology faster and more profitable.

Aesthetic Press has developed a groundbreaking change in how porcelain systems are designed.

The AP press ceramic (feldspar based) can be re-pressed onto an already pressed crown or bridge. The option of layering and the advantage of re-pressing ingots, which have already been used, is just another example of how AP is both a practical and economic system.

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
Silicate 400 MPa

General Overview

The APVolution M porcelain offers a wide range of options for technicians using a layering or press over technique onto a metal frame work.

Technicians can chose between a very economic, streamlined and production oriented selection of powder porcelains and ingots as well as an assortment for the high end technician looking for a broad spectrum of modifiers and special effects powders.

The assortment offers High Fusing as well as Low Fusing Powders, which compliment each other and offer therefore a perfect solution for your daily cases at the bench.

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PRESS OVER METAL INGOTS

Press & Stain / Press & Maltechnik

Save money on inventory by applying the art of stain techniques to one of our value based ingots (light, medium, dark, or dark plus) and corresponding powders. Use one of our four chroma shades (page 48) along with effect stains to create the perfect match for your patient. One of the significant advantages is that the application of glaze paste is not necessary. A natural glaze ensures a homogenous and highly aesthetic result.

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Press & Layer / Press & Schichttechnik

The Aesthetic-Press Dentin ingots are mainly used for anterior restorations. After transferring the wax up to a full contour body shape, a precise cut back defines the dentin core. The Aesthetic-Press ingots are designed to match the desired chroma hue and translucency together with the AP Stain Kit. There are 14 enamel and translucent powders available to complete the restoration to its desired form, shape and color effect.

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PRESS OVER METAL INGOTS

Gingiva

There are no limits for pressing extensive reconstructions with "pink" and "white". Even 14 units can be made with the Aesthetic-Press ceramic without any problems. The Double Press -Technique is explained on our website www.apdental.net.

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Layering / Schichttechnik

These build up powders can be conveniently used to either layer or to correct any pressed crowns. The blended powders exactly match the colors of the ingots. Therefore, any repair or added powder will blend in nicely for an aesthetically pleasing result.

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We have now expanded our line to include Low Fusing Powders for porcelain to metal. This can be used if a porcelain margin has been pressed to maintain a sharp and crisp margin.

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10

APVOLUTION M

For the value based ingots, we recommend the following opaque shades:
 Opaque Light for A1, A2, B1, B2, C1, C2, D1, D2
 Opaque Dark for A3, A3.5, A4, B3, B4, C3, C4, D3, D4
 Opaque Bleach for Bleach Shades

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Working with AP-Efficiency Through Simplicity

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General Overview

The APVolution Z porcelain, offers a wide range of options for technicians using a layering or press over technique onto a zirconia framework.


Technicians can choose between a very economic, streamlined and production oriented selection of powders porcelains and ingots as well as an assortment for the high end technician looking for a broad spectrum of modifiers and special effects powders.

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PRESS OVER ZIRCONIA INGOTS AND POWDERS

Press & Stain / Press & Maltechnik

Save money on inventory by applying the art of stain techniques to one of our value based ingots (light, medium, dark, or dark plus) and corresponding powders. Use one of our four chroma shades (page 48) along with effect stains to create the perfect match for your patient.

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APVolution Z Value Based Ingots

**Please see page 33 for layering powders

Product Name/ Artikel	Product Code/ Artikelnummer	Product Code/ Artikelnummer
Size/Größe	2g Ingots/Pellets	3g Ingots/Pellets
APVOLUTION Z LIGHT	2ZL	3ZL
APVOLUTION Z MEDIUM	2ZM	3ZM
APVOLUTION Z DARK	2ZD	3ZD
APVOLUTION Z DARK PLUS	2ZD+	3ZD+

PRESS OVER ZIRCONIA INGOTS

Press & Layer / Press & Schichttechnik

APVolution Z Dentin Ingots are designed for anterior and posterior restorations for specific shades with minimal staining. These are available in 16 main VITA color categories: A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4 and three bleach shades. With the dentin ingots, the technician can do a cutback, internal staining and add enamel and transpa layers. 14 Enamel, Opal and Transpa Powders are available.

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General Overview

The APVolution S ingots is a highly durable, pressable lithium silicate porcelain with flexural strength of 400 MPa +- 40.

In combination with the newly patented APColor Key Shade Guide, this product is unique and offers significant advantages over regular shade guides.
See article.*

A wide range of shades and opacities are available to match even the most difficult shades.


The APVolution S ingots are available as LT, MT, HT, HT+ and Bleach ingots.
The new APColor Key Shade Guide is made out of the same ingots, which the technician can now use for their next case.

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400±60 MPa Lithium Disilicate

Color Key Shade Guide

The only shade guide
on the market with
12 stump shades

Indicator + Stump = Shade

Realistic indicator thickness of 0.9-1.1mm

LT, MT, HT, HT+ & Bleach Shades

Find the right opacity and value before you start the case

www.apdental.net

Product Name/ Artikel	Product Code/ Artikelnummer
Size/Größe	3g
APVOLUTION S LT A1	SLTA1
APVOLUTION S LT A2	SLTA2
APVOLUTION S LT A3	SLTA3
APVOLUTION S LT A3.5	SLTA3.5
APVOLUTION S LT A4	SLTA4
APVOLUTION S LT B1	SLTB1
APVOLUTION S LT B2	SLTB2
APVOLUTION S LT B3	SLTB3
APVOLUTION S LT B4	SLTB4
APVOLUTION S LT C1	SLTC1
APVOLUTION S LT C2	SLTC2
APVOLUTION S LT C3	SLTC3
APVOLUTION S LT C4	SLTC4
APVOLUTION S LT D2	SLTD2
APVOLUTION S LT D3	SLTD3
APVOLUTION S LT D4	SLTD4

Indication: Full crowns and Veneers. In order to get a natural looking crown, the MT ingots can be used for both the press/cut back /layering technique and press & stain technique.

The technician can define the amount of opacity needed through a defined depth cut - internal staining and micro-layering.

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APVOLUTION S INGOTS

LT - Low Translucent

The LT ingots show a high level of chroma and opacity, which helps to cover dark stump shades or certain discolorations.

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Indication: Full crowns & bridges. In order to get a natural looking crown, the LT ingots can be used for the press/cut back and layering technique. The additional transparent layer will give the crown a life-like appearance, yet the core material will ensure the chroma needed to match the shade.

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APVOLUTION S INGOTS

MT - Medium Translucent

The MT ingots show a medium level of chroma and opacity, which helps to cover light stump shade discolorations.

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The HT ingots show a high level of translucency.

Indication: In order to get a natural looking crown, the HT ingots can be used for the press and stain technique. The overall appearance will be a natural translucent looking crown or veneer or even for non prep cases.

The HT+ ingots are designed for thin veneers or inlay only restorations offering a maximum of translucency.

Product Name/ Artikel	Product Code/ Artikelnummer
Size/Größe	3g
APVOLUTION S HT + 1	SHT+1
APVOLUTION S HT + 2	SHT+2
APVOLUTION S HT + 3	SHT+3

Product Name/ Artikel	Product Code/ Artikelnummer
Size/Größe	3g
APVOLUTION S HT A1	SHTA1
APVOLUTION S HT A2	SHTA2
APVOLUTION S HT A3	SHTA3
APVOLUTION S HT A3.5	SHTA3.5
APVOLUTION S HT B1	SHTB1
APVOLUTION S HT B2	SHTB2
APVOLUTION S HT C1	SHTC1
APVOLUTION S HT C2	SHTC2
APVOLUTION S HT C3	SHTC3
APVOLUTION S HT D2	SHTD2
APVOLUTION S HT D3	SHTD3

APVOLUTION S INGOTS

Bleach Ingots

The eleven bleach shades are available in different opacities as LT, MT and HT ingots offer a wide spectrum to choose from.

Indication: In order to get a natural looking crown, the Bleach ingots can be used in a variety of scenarios for either press & stain, or cut back & layering technique. The technician can even increase the value if needed by adding a thin layer of transpa white layering powder mixed with a white modifier. The APVolution S&Z Powder selection offers a wide range of powders compatible with the APVolution S ingots.

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
Temperatures can vary between different furnaces up to 40°C.

Observe the visual result and contact customer support in case of questions.

FIRING CHART - BUILD UP POWDER / BRENNTABELLE - SCHICHTKERAMIK

S & Z
Powder / Pulver

Build Up Powder	Idle Temp	Rising °C/min	End Temp	Preheat	Holding Time	Vac on
Schichtkeramik	Bereitschaft	Steigrate °C/min	Endtemp.	Vortrocknen	Haltezeit	Vak an
Dentin / Incisals	400°C	45°C	780°C	4 min	1 min	450°C
Glaze Bake / Glanzbrand	400°C	45°C	750°C	4 min	1 min	-
Correction Powder / Korrektur Pulver	400°C	45°C	720°C	4 min	1 min	450°C

APCOLOR KEY SHADE GUIDE

SHADE GUIDE

APColor Key Shade Guide

Product Name/ Artikel	Product Code/ Artikelnummer
APVOLUTION S LT	APCKLT
APVOLUTION S MT	APCKMT
APVOLUTION S HT & HT +	APCKHT
APVOLUTION S BLEACH	APCKBL
APVOLUTION S STUMP	APCKST

APVolution S Color Key Shade Guide

Unlocking the mystery of shade communication

By Jörg Müller, MDT, CEO of Aesthetic-Press

INTRODUCTION

Taking the right shade for a patient is the everyday challenge for dentists and for technicians. To find the right matching shade, many variables have a big impact on the result. Some of them are the light source, the time of day, the angle at which the light shines onto the tooth, the surface texture and reflection of the tooth, the density of the material and the underlying stump shade or framework for the final restoration. The standard shade guides offered by the industry up to now, have certain unrealistic features technicians have difficulties to match.

First, the thickness of the shade guide is sometimes 3 to 4 x thicker than the actual crown or veneer restoration. Since color is a result of various layers, shade matching depends on the skillset and artistic work of the technician to recreate the effects, shade, value and chroma. In simple words, the technician has to work his magic to recreate a color in a crown which oftentimes is only 0.6 mm - 1.2 mm thick vs. the shade guide showing a thickness of up to 4 mm.

THE CURRENT PROBLEM

Shade guides are often:

- too thick
- porous surface
- made out of a material the patient case is not going to be made out of
- impossible evaluation of opacity level
- impossible evaluation of stump shade effect to the desired final shade.

With the current porcelain systems on the market, the technician is not able to see or evaluate how the color of the porcelain restoration will be influenced by the underlying stump shade before starting the case. The stump shade is the most critical aspect in any full porcelain restoration, either a veneer or a full porcelain crown. The stump shade influences the final color and often times has a negative effect on the shade of the crown.

This can only be avoided by choosing a certain material opacity or value to block a dark or grey looking stump. Therefore, most of the shade guides of full porcelain systems on the market offer color indicators which show a porous and unfinished surface. This leads to an unsatisfying evaluation of the true shade and most importantly does not allow different case scenarios with the variety of different stump shades.

SOLUTION

...see the shade before you start the case!!

The design of the AP Color Key Shade guides shows a thickness of 0.9 mm at the incisal edge, about 1.00 mm in the center of the indicator and around 1.3 mm towards the cervical or neck of the tooth. These thicknesses are realistic for a technician to match a shade, as they reflect the average thickness of a crown.

WHY SHOULD DENTISTS USE THE AP COLOR KEY SHADE GUIDE?

The current shades guides on the market only offer a 2 dimensional approach to color and show only chroma and value. In terms of opacity, they are only 1 dimensional.

Shade guides show a color, for example a VITA A2 only in 1 level of opacity, whether the patient is 18 or 80 years old. The underlying shade can be an A2, yet the depth and the translucency are completely different, an effect the generic shade guides are not able to match.

The APColor Key Shade Guide offers a 4 dimensional approach.

- Color in A, B, C and D Chroma,
 - Value in 1, 2, 3 and 4
 - Opacity, eg LT, MT, HT and HT+ in combination with the stump shade
 - Stump shade in 1-8
- Choose a shade with the material the technician is working with.
 - Choose a shade with a realistic thickness of indicator.
 - Choose the right kind of indicator in conjunction with the stump shade that comes closest to the patients shade.

ADVANTAGES

- Higher success rate-Less re-do rate. Select matching ingots before the case begins
- Traditional shade guides mostly only offer shade selection in one opacity, missing a big part of the color issue. Choose matching opacity/translucency: MT, LT, HT, HT +
- Shade guide manufactured with the porcelain available for the case - ensures a more realistic color selection
- Technicians can select ingots with their eye rather than depending on an app
- Much improved communication between dentist and patient, chairside, true color selection possible.

APColor Key Shade Guide Indicators

Process in pictures

SHADE GUIDE SELECTION

The APVolution S Color Key Shade Guide offers 4 different translucencies/opacities. The technician can choose which ingot is most suitable to match the patient's shade according to their stump shade.

1. Low Translucency or LT

This material shows more opaque porcelain which can cover dark stump shades or certain discolorations.


Indication: In order to get a natural looking crown, the LT ingots can be used for the press/cut back and layering technique. The additional transparent layer will give the crown a life-like appearance, yet the core material will ensure the chroma needed to match the shade.

2. Medium Translucency or MT

This material shows a medium opacity level to cover stump shades with mild discolorations.

Indication: In order to get a natural looking crown, the MT ingots can be used for the press/cut back and layer technique, or press and stain technique. The technique applied will depend on the technician's visual evaluation.

3. High Translucency or HT and HT+

This material shows a high translucency level to cover stump shades with no discoloration.

Indication: In order to get a natural looking crown, the HT ingots can be used for the press and stain technique. The overall appearance will be a natural translucent looking crown.

The HT+ ingots are designed for thin veneers or inlay only restorations offering the maximum translucency.

4. Bleach Shades

The eleven bleach shades are available in different opacities such as LT, MT and HT ingots and offer a wide spectrum to choose from.

Indication: In order to get a natural looking crown, the Bleach ingots can be used in a variety of scenarios for either press & stain, or cut back & layer technique.

The technician can even increase the value if needed by adding a thin layer of transpa white layering powder mixed with a white modifier.

The APVolution S&Z Powder selection offers a wide range of powders compatible with the APVolution S ingots.

5. Stump Shades

The stump shades are made in a shape to match the concavity of the APColor Key Shade Guide indicators. To get the best possible light transmission, a thin film of Vaseline can be applied onto the Stump Shade indicator to ensure the best possible light transmission into the porcelain indicator.

...UNLOCK THE
MYSTERY OF SHADE
COMMUNICATION

INNOVATION:

Innovation: The concavity of the manually finished and polished porcelain indicator in combination with the stump shades offer an excellent way to select the right porcelain shade and opacity before starting the case.

The perfectly finished surface texture of the indicator enhances the reflection and visual result comparable to a high end restoration.

This increases the efficiency and success rate of any full porcelain restoration.

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...UNLOCK THE MYSTERY OF
SHADE COMMUNICATION

APVOLUTION **S**

Silicate 400 MPa

APVOLUTION **Z**

Press to Zirconia

LAYERING AND CORRECTION POWDERS

YOU CAN LAYER ON BOTH!

S & Z Powders / S & Z Pulver

The APVolution S & Z powder can be layered onto zirconia frames and silicate / lithium silicate crowns

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APVOLUTION **S & Z** POWDER / PULVER


S & Z Gingiva Powders / S & Z Gingiva Pulver

APVolution S & Z gingiva powders are offered in 8 different shades to match the natural gingiva. These powders can be applied to any zirconia frame and are compatible with the APVolution S & Z ingots and powders.

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There are 16 dentin powders, 3 transpa dentin, 14 enamel, transpa and opalescent powders available as well.

The dentin powders can be used for regular build up as well as for correction purposes.

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



Where Form and Function Meet

Using the Aesthetic-Press system to create a full-arch zirconia bridge.

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Over the years, I have seen many restorations being modified and compromised due to poor planning from the start. Developing a plan, turning it into a protocol and then following every step precisely, will ensure a predictable and successful restoration.

For aesthetic anterior restorations it is most important to start with a diagnostic wax up. This is a solid foundation for any case. All the critical points such as length, proportions and anatomical details will be determined. No technician wants to make changes – such as adding to the incisal edge – when it's time to deliver a case. We feel great pleasure, when the case just drops in and we meet the patient's expectations. To make this a standard rather than playing a lottery game, we need to follow a step-by-step approach.

First Step – the Wax Up: Find the right functional elements, the right size, proportions, and horizontal and vertical dimensions, which are parameters that can be determined in the wax up.

OCCLUSAL CONCEPTS

For edentulous patients, it has been common practice to set up denture teeth because they are easy to place and quickly lead to a result. However, because those denture teeth are more generic rather than customized, it may lead to limits if the form or function needs to be modified.

A more suitable solution, and equally as efficient, is the OccluMaster wax patterns by Aesthetic-Press. Because of their flexibility, these fully contoured posterior crowns help find the right occlusal dimensions in a very short time, while adjustments or functional elements can be easily done. With the OccluMaster, a variety of occlusal concepts can be utilized and easily modified. The conventional P.K. Thomas anatomy can be found in our classic design, and the more advanced theories, like the Biomechanics taught by the renowned Michael Polz from Germany, can be found in the Premium morphology (Fig. A).

CREATING THE RESTORATION

Milling a bar in zirconia offers a variety of advantages compared to its counterpart, metal. The very minimal porosity level of zirconia is far more biocompatible than any metal construction, not to mention the implications of soldering connections.

A passive fit of the bar is essential, which will be verified through a try in, in its frame stage. Any adjustments can be done chair-side with light cure materials. With a thorough milling procedure and controlled sintering process, the bar fits nicely to the implant.

Fig. A OccluMaster Premium wax patterns can help set the occlusion at the start of the case.

Fig. B The framework is milled from the zirconia block.


Fig. C The framework is colored and fired while still connected to a piece of the block.

Fig. D The polished zirconia bar seated on the model.

Fig. E Light cure material is applied to the frame before it is placed in the matrix created from the original wax up.

Fig. F For the "drop the pin" milling technique, the milling bur is dropped approximately 0.75mm...

01 The bar is designed and milling is started using a 16mm zirconia block. Four different burs are used to mill the bar on a Zirkozahn Eco milling machine (zirkozahn.com). The abutment arias need to be milled very carefully.

02 Once the bar is almost finished, the housing for the screws is drilled into the zirconia before sintering. The Eco milling machine provides the ability to mill the housing for the lateral screws, especially for the constructions like this complex restoration (Fig. B). In addition, the bar can be milled either parallel or at a 2° angle. This will ensure that the over structure will fit nicely, and the friction will keep both parts together. The set screws will lock the construction in its final position. The 4 lateral screws carry the metal housing, which will be inserted once sintering is completed.

03 After milling the bar, the color liquids are added and the bar is sintered while still connected to the block to retain the precise shape during the sintering process (Fig. C).

04 Cutting the connectors had to be done very carefully with water-cooling to avoid a fracturing of the bar. After fitting the bar, the zirconia was polished with Oscar, AP Diamond Polishing Paste, which has a high diamond content and gives the bar a perfect surface shine (Fig. D).

05 To cover the bar of the overstructure, the Rigid light cure material was applied before the frame was put into the silicon matrix which was created from the original wax up (Fig. E).

06 Maintaining the occlusion is a very critical point when duplicating a precise wax up. A special mold maker is ideal for transferring the wax up into resin. This mold maker can also be used for a wax injector to re-inject the previous wax up.

07 After milling the over structure, the final details were milled after utilizing the “drop the pin” technique. The milling bur is dropped approximately 0.75mm, whereas the tester will stay in its original position (Fig.F). This technique will ensure the pressed porcelain has enough room for maximum stability and perfect light reflection (Fig. G).


Fig. G ... ensuring the milled zirconia leaves appropriate room for the pressed porcelain.

08 Using the translucent Aesthetic-Press Zircon ingot will give you the optimum porcelain support. As Figure H illustrates, the Eco milling machine can be set to different angles in order to mill every corner of the restoration. After opening the set screws of the hinge, the milling plate can be adjusted as needed. The rotation of the 5 axes allows the placement of the lateral screws.

09 The Gingiva, as well as the tooth colors, are applied carefully with a brush before the sintering process. The color liquids are needed to support the final color with the right Chroma from within. For the experts, the slight color nuances can be applied with a micro brush. A variation from violet, grey as well as a mild orange can be applied to bring some life into the zirconia structure (Fig. I).

Fig. H The Eco milling machine can be set to different angles to mill every corner of the restoration.

10 The pink porcelain will be pressed first. To accomplish this, the original wax up is re-injected onto the framework. The AP injection wax, with its specific properties, allows the occlusion and crowns to be removed with a simple cut of a scalpel at the red/white junction allowing the wax to be peeled off.

11 The length of the sprues for pressing this porcelain were significantly longer than allowed in any textbook. In order for a successful press result, the restoration was supported by 4 sprues on the inside and 5 on the outside of this full arch zirconia framework (Figs J, K).

12 After investing the bridge for 15 minutes, bench set time was used before placing the ring into a cold burnout furnace. The program was set with a rising temperature of 17 degrees Celsius per minute. After reaching the final temperature of 850 degrees Celsius, a holding time of 90 min needs to be completed before placing the ring into the press furnace. Pressing a large restoration as described, the press temperature should be increased by 40° C to the general temperature listed in the firing chart.

13 The AP Feldspar ingots should not be preheated. Using 6 x 2.5g ingots, the length of the plunger had to be shortened by 50% in order for the oven to close. The gingiva part was pressed with the Aesthetic-Press gingiva ingots. The press time was 20 min and the distance travelled was 21mm. The press time describes the time used for the plunger to move the porcelain into the form. Working with such a large amount of zirconia requires a lot of attention to a very slow rising temperature and cooling time. An extended cooling time was needed before carefully de-vesting the massive zirconia overstructure. After checking the fit of the upper construction with the model, the wax injection completed the missing layer of the clinical crown.

Fig. I Gingival and tooth colors are applied to the framework to support the layered porcelain.

Fig. J & K Extra long sprues are used to properly press the pink porcelain layer.

Fig. L Each tooth is connected by its own sprue when the white porcelain layer is pressed.

Fig. M The complete and glazed restoration.

Figs. N & O The wax up is duplicated in the final result with the planned occlusion precisely fabricated.

Fig. P The complete restoration in two pieces.

14 Little adjustments were needed before the construction is ready to be re-invested for the second press cycle to press the white porcelain. Each tooth was connected with one sprue of gage 8 (**Fig. L**). The reconstruction should be centered as much as possible to keep a balance of the porcelain during the cooling phase. This proves the point, that with AP ingots, the structures of any kind and size can be successfully pressed. The press temperature and thickness are the keys to success.

15 With a precise wax up, only minimal corrections are needed to complete and glaze the construction (**Fig. M**). The ability to finalize the occlusion in wax makes such a large and complex case much easier to plan and manufacture. All the functional elements can be designed at the right place, more easily and more precisely than with a conventional build up technique in porcelain (**Figs. N, O**).

CONCLUSION

At the time when I developed the Aesthetic-Press System, I was looking for a system, which would help every technician simplify their work and make it more predictable at the same time (**Fig. P**). Working in my laboratory with my clients and technicians, enabled me to develop this proven system and make it a success on a daily basis. This system is a great opportunity for every technician to bring high-end quality at an affordable price to their laboratory and even more importantly to make your work consistent for you and your technicians case by case.

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
"So good it deserves an Oscar!"
- Jürgen Auffarth, Salzburg

POLISHING PASTE / POLIERPASTE

Use this polishing paste for all your polishing needs. Ensure a quick, homogenous result without smearing.
 "This polishing paste deserves an Oscar!" - Jürgen Auffarth

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The ideal polishing paste for a perfect high shine!

The high diamond content of the "Oscar" Diamond Polishing Paste ensures a quick result for all surfaces without smearing. Composite, metal, ceramic and even zirconia surfaces can be treated with this phenomenal product.

Every technician will enjoy the ease of use and quickness of its result.

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EFFECT STAINS / EFFEKT MALFARBEN

Effect Stains / Effekt-Malfarben

These stains are low fusing and fluorescent paste stains. This paste material can be applied in very thin layers to achieve the most difficult color effects.

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For pressable ceramics, paste stains are recommended. Depending on the effect, 1-2 glaze cycles are recommended. These Low Fusing paste stain enables a nice surface texture and shine at 760 Celsius (1400 Fahrenheit).

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With the Chroma Stain shades A,B,C,D the desired chroma shade can be achieved for the "Classic" or the "Zircon" ingot.

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MARGIN STAIN

The Margin Stain is one of the most important colors within the stain kit. With the more opacous color, transparent porcelain margins can be covered by using the press on metal technique.

For anterior restorations, the Margin Stain can very nicely mimic mamellon effects.

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GLAZE PASTE / GLASUR-MASSE

The glaze paste is a low fusing material used to glaze veneers or porcelain inlays. This glaze paste leads to great results applied in thin layers, a preheat temperature of 450 Celsius (842 Fahrenheit) and a 6 min dry time. At times, two glaze cycles are recommended.

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WAX INJECTOR / WACHSINJEKTOR

This wax injector is very user friendly. The air pressure inside the injector is built up using a simple hand pump. The pressure should be between 0.5 and 1 bar. The temperature should be set around 74 Celsius.

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With this easy to use technique, even large restorations can be duplicated and reproduced very easily.

We can achieve a 1:1 translation of the wax up to the definite restoration, which makes every case plan and diagnostic wax up a successful case strategy!

This technique saves a lot of time and effort to reproduce the diagnostic wax up done in the planning stage of a case. With this procedure we demonstrate efficiency and accuracy while guaranteeing a controlled outcome with each case.

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
AP Silicone Blue Standard Pack Prod. Nr. 7100
base, jar of 800 g
catalyst, jar of 800 g
2 measuring scoops

Mega Pack Prod. Nr. 7101
base, bucket of 5 kg
catalyst, bucket of 5 kg
2 measuring scoops

AP Silicone Blue – 86 Shore A
very high consistency, easily kneadable, precise detail reproduction, stable in shape and volume

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AP Silicone Blue Rock Standard Pack Prod. Nr. 7102
base, jar of 800 g
catalyst, jar of 800 g
2 measuring scoops

AP Silicone Blue Rock – 96 Shore A
high viscosity, smooth initial consistency, high definition, dimensionally stable

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WAX FLAKES / WACHS FLAKES

This specially designed injector wax is very suitable for all pressable crowns and bridges. Its ash free content, burns out nicely to ensure beautiful vital results.

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Workflow: After establishing a diagnostic wax up and with the help of a silicon matrix, the wax injector will help reproduce this exact situation by injecting the wax onto any framework or prepped model.

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The composition of this injector wax allows very thin areas to flow nicely and its carving capabilities are extremely easy.

Use this wax to make perfect adjustments every time to any form you are working with.

A diagnostic wax up can easily be transferred to the definite case using a matrix. By doing this, the technician saves a lot of time avoiding a re-wax. Also, any temporary that has been established in the patient's mouth, which reflects all aesthetic and functional requirements, can now be translated without any compromise into a beautiful long lasting restoration.

The beige color is very pleasant to the eyes and the added vanilla scent makes this material a pleasure to work with. This wax is available in 500g.

Product Name/ Artikel	Product Code/ Artikelnummer
Size/Größe	500g
WAX FLAKES / INJECTOR WAX	9558

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OCCLUMASTER KITS

With the OccluMaster, we have created a new dimension in dental restoration, ensuring a perfectly stabilized and balanced occlusion.

The properties of the OccluMaster are ideal for handling large and even small reconstructions. Available in 4 different styles.

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OccluMaster Kit Contents

16 Sizes (10 each) / 10 mal 16

This wax can be carved like regular carving wax and contours can be modified very easily. By adding warm air from a hairdryer, the wax becomes elastic and works well for cases with even little reduction.

We combined all the functional elements that lead to a physiological restoration according to the Bioaesthetic and Biomechanical principle.

Since the OccluMaster is a standardized form for each posterior tooth, each and every tooth will be the same and guarantee a maximum in consistency.

With one OccluMaster you buy more than two decades of accumulated knowledge in form and function.

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OCCLUMASTER REFILLS

OccluMaster Classic Refills

The wax properties allow the technician to place the prefabricated full anatomic wax crown into occlusion without breaking of the wax pattern.

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

Working Time

Using the OccluMaster	
Wax up	10 min
Invest and deinvest per unit	3 min
Fit check, finish & glaze	15 min
Working time on porcelain margin:	0 min
Total	28min

...Compare

Working Time

Regular Build Up	
Working time on build up technique	50 min
Working time on porcelain margin 2 additional firings plus finishing and fitting:	20 min
Total	70 min

Gain of 150% in time
from start to finish

OCCLUMASTER PATTERN GUIDE

OccluMaster Classic

OCCLUMASTER PATTERN GUIDE

OccluMaster Premium

OCCLUMASTER PATTERN GUIDE

OccluMaster Natural

OCCLUMASTER PATTERN GUIDE

OccluMaster Classic Plus

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Workflow: After establishing a diagnostic wax up and with the help of a silicon matrix, the wax injector will help reproduce this exact situation by injecting the wax onto any framework or prepped model

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8 upper & 2 lower forms (no wax patterns)

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Anterior Master Kit 2 Forms + Patterns

2 forms of choice and the corresponding wax patterns (and example is shown below)

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







Using the Aesthetic Press Vanilla Wax in a multi-quadrant case.

by Jörg Müller, MDT


For dental technicians, the move to all digital has not always been very easy. The skill to work with a model, with wax or porcelain, or any number of manual skills that were earned over a long period of time are certainly more than planning and designing at a computer screen. Trying to see and work in three dimensions even though there are only two dimensions on a screen can be a challenge.

PLANNING DEPARTMENT

A thoroughly planned case always starts with a wax up. Comprehensive cases such as multiple quadrants or a full mouth reconstruction need precise planning. But we don't even have to go that far. Any anterior restoration with more than two units at least should be waxed up in order to find the right length and position of the teeth. This wax up should be translated into a temporary for the patient to review and get comfortable with the new esthetic smile designed by the team effort of the technician and dentist. We can achieve this workflow either manually or digitally. However, not every software library produces teeth in a "nature"-like design. Oftentimes, the milled units coming out of the milling machine need some touch-ups. Either it is the marginal area that might need some attention, using a coping for casting metal, or using a full contoured crown for press a porcelain crown. Next to the marginal areas, details like occlusal anatomy or functional element might need some reshaping.

A MULTI-PURPOSE WAX

Most wax blanks have a higher content of resin filled particles, which means the "wax puck" will be very rigid and harder to make any adjustments, not to mention a sealed margin. A crown, which is milled out of a material that handles like conventional carving is ideal to modify. Aesthetic Press Vanilla Wax is specially designed to serve multiple purposes. One can use it as wax flakes for a wax injector, as regular carving wax and also as a wax blank for most CAD/CAM systems on the market. This level of multitasking in a wax is unique in the industry. Another aspect which is pleasant for every technician—besides the vanilla scent, of course!—the color of the wax is easy on the eyes and as a technician I can finally "see the anatomy." This wax is even able to be used for an esthetic try in, or for functional control movement for posterior restorations.

PRACTICAL CASE

1. After scanning and designing the multi quadrant case digitally, the case is milled out of the 98mm AP Vanilla wax blank.
2. The anatomy is fine; however, certain details are missing. These details are not always easy to design at a computer screen and the given tools of the current software available. With the help of an appropriate carving instrument—e.g. the PTC Wax Carver No. 1—the anatomy can be refined to its desired functional design.
3. Certain functional elements are more easily checked and designed in a physical articulator. The split file technique is today's frequently used workflow to effectively and precisely plan and fabricate a porcelain restoration.
4. The right framework design is essential for the durability of any restoration. Designing the full contour crown and being able to reduce the design by approximately 0.8 mm serves enough space for a nice esthetic result but also for a stable occlusion, minimizing the risk of chipping. Press to zirconia
5. After pressing the units with the AP Zircon for Staining ingots, very little refinements are necessary.
6. The feldspathic porcelain allows easy repair or adjustment if necessary. The main advantage here is that the glaze cycle is done without the use of a glaze paste. The self-glazing properties save most technicians a lot of time to stain the right shades to the porcelain crowns.

CONCLUSION

In my experience, combining the best of both worlds is a perfect connection of the dots. Being able to manually modify a CAD CAM supported restoration with a material that enables the properties a technician is used to maintains the joy and fun factor at work. A seamless integration of material and workflow is the most important guarantee and success for each case, but also maintains the happiness of each technician!


Figure 1

Figure 2

Figure 3

Figure 4

Figure 5

Figure 6

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VANILLA WAX BLANKS

98mm for Open Systems

This true wax blank is the ideal material to mill wax. With its unique properties the milled crowns can be easily adjusted.

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110mm for Dental Mill Machines

**Important
Information**
**Wichtige
Information**

Important Parameters:
Drill Speed – 6000-9000 RPM
Feed Rate - 9mm/sec
Burs – PMMA or Zirconia burs no more than 2 flutes

TIP! Refrigerate blanks for extra firmness.

Size/Größe	Product Code/ Artikelnummer
110 x 12mm	9567
110 x 14mm	9559
110 x 16mm	9560
110 x 18mm	9561
110 x 20mm	9562
110 x 25mm	9563
110 x 30mm	9564

VANILLA WAX BLANKS

Wax Blank Milled for Jet Blank

Size/Größe	Product Code/ Artikelnummer
Milled for Jet Blank	9565

AP VANILLA SPEED WAX BLANKS

Vanilla SPEED 98mm for Open Systems (NEW FORMULA)

This new formula of wax has easy modification similar to our original Vanilla Blank with fast, accurate milling capability. It's the perfect blend of wax and resin material with the same pleasant scent.

For a fast mill you can work with, try our new Vanilla SPEED!

Size/Größe	Product Code/ Artikelnummer
APVanilla SPEED 98 x 12mm	9580
APVanilla SPEED 98 x 14mm	9581
APVanilla SPEED 98 x 16mm	9582
APVanilla SPEED 98 x 18mm	9583
APVanilla SPEED 98 x 20mm	9584
APVanilla SPEED 98 x 25mm	9585
APVanilla SPEED 98 x 30mm	9586

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AP VANILLA SPEED WAX BLANKS

White Lemon 98mm for Open Systems (For Diagnostic)

A beautiful white wax blank perfect for diagnostic wax ups. Fresh lemon scent and easy to work with.

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APCAD Blue Blocks- Lithium Disilicate milling blocks. The high strength performance and amazing aesthetics of this material make it our no.1 choice for restorative dentistry.

The APVolution Color Key Shade Guide makes this product unique:

- Choose a shade with the material the technician is working with
- Choose a shade with a realistic indicator thickness
- Choose the right kind of indicator in conjunction with the stump shade that comes closest to the patients shade
- modify the indicator while taking the shade with your stains

Shade guide available in 1mm and 0.5 mm thickness for crowns, veneers and minimal invasive cases

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FIRING CHART

FIRING CHART / BRENNTABELLE

CRYSTALLIZATION PROGRAM - APCAD S MILLING BLOCKS

	Idle Temp Bereitschaft	Dry Time Vortrocken	Rising °C/min Steigrate °C/min	End Temp Endtemperatur	Holding Time Haltezeit	Vac Start	Vac End
Phase 1	400	6 min	90° C	820° C	0:10	550	820
Phase 2	–	–	30° C	840° C	7:00	550	840

Easy-Vest is a high quality, high-speed investment for crown and bridge work, as well as for all pressable systems.

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The unique ability to use this material for all purposes enables laboratories to reduce inventory and overhead.

Easy-Vest is an all-in-one carbon free phosphate bonded investment:

- Casting alloys
- Pressable porcelains
- Suitable for all types of dental alloys and pressable porcelains
- Ringless techniques for high speed and slow speed techniques possible
- Optimal expansion can be controlled with the Easy Vest Liquid
- Compatible with non precious and Implant substructures

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EASY VEST PREMIUM

The Easy Vest Speed Premium is primarily made for press ceramics to offer a perfectly smooth surface without any reaction layer especially for Lithium Disilicate materials. Fine grain, easy to divest.

Product Name/ Artikel	Product Code/ Artikelnummer
EASY VEST PREMIUM 5kg INCLUDING LIQUID	9010P
EASY VEST PREMIUM LIQUID 1 LITER	9012P
EASY VEST PREMIUM POWDER ONLY	9011P

No Reaction Layer!

Specifically for Lithium Disilicate

PLUNGER MAKER / PRESS STEMPEL FORMER

With the tear free silicon material, individual press plungers can be poured with the leftover Easy Vest investment. Using the AP plunger maker will lower the cost of investment by 50%.

Each time a ring is invested it can produce up to 3 plungers. After 60 min the excess material can be trimmed with the model timmer. Push the plungers out from the bottom and use compressed air to blow them out of the form. Preheat plungers to eliminate excess moisture.

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
...MAKE YOUR OWN PLUNGERS AND
SAVE \$\$


INVESTMENT RING / MUFFEL

The AP 200, 300 and 400 are designed with many special features. The tapered form makes it easy to release the investment out of the rubber mold, which reduces the risk of micro fractures that lead to ring failure.

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The AP 200, 300 and now the new 400 are the only investment rings on the market with a number system on the inside. The numbers show the exact position of the crowns after the press cycle and enables a quick de-vesting process. The lines show precisely, where to cut the ring in order to break off the excess investment parts. The numbers can be found on the bottom of the molds as well as the inside of the silicon ring and on the inside of the lid. The ring can only be closed (or put together) in a precise position due to the individual locks. The piston of the AP 200 is built higher than conventional rings, which allows pressing 3, 2g ingots.

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AQUALINE SMART

In a compact design, the Aqualine Smart offers everything a perfect ceramic workplace needs. The moisture of the materials is constantly kept at the perfect level via cellulose strips, without producing any unpleasant odors.

There is no chance of materials drying out. It saves on cost and materials. The overall ergonomic design allows efficient working processes in the smallest of spaces - simply smart, isn't it?

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Aqualine Smart tray
Color: Anthracite rubber touch
Dimension: 260 x 170 x 32 mm
Components: tray with glass plate and hinged cover, self adhesive write on labels, cellulose strips

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AQUALINE SMART POOL SET

Aqualine Pool Set
Dimension:
105 x 65 x 48 mm
Components:
250 ml rinsing cup, sponge and sponge holder

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AQUALINE SMART STAIN TRAYS

Aqualine Stain Tray
Color:
White porcelain
2 Piece Design:
1 stain tray, 1 mixing tray

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The Aqualine Mini is specially developed for smaller working and ceramic repairs. The very compact tray subdivided into three segments enables dental technicians to perform efficient layering and keep ceramics perfectly moist. Fits in the smallest of places. (13.3 x 7.1 x 2.3 cm)

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AQUALINE MINI STAIN TRAY

Aqualine Mini Stain Tray provides professional ceramists with two modern mixing bases in miniature format, which constantly guarantees the best modeling and material properties as well as brilliant firing results - without the need to constantly rewet and spatulate the mixed ceramics and without the risk of drying out or contamination.

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Why customers chose US!

- 20 years experience
- Competitive prices
- Unmatched service

AUSTROMAT 624 by DEKEMA

AUSTROMAT 654 press-i-dent by DEKEMA

Product Name/ Artikel

Product Code/ Artikelnummer

AUSTROMAT 624 by DEKEMA	9134
AUSTROMAT 654 press-i-dent by DEKEMA	9102
AUSTROMAT 664 by DEKEMA	9135
AUSTROMAT 674 by DEKEMA	9136
AUSTROMAT μ SIC by DEKEMA	9137
AUSTROMAT 624 i by DEKEMA	9130
AUSTROMAT 664 i by DEKEMA	9132
AUSTROMAT 674 i by DEKEMA	9133

See website for more information

AUSTROMAT 664 by DEKEMA


AUSTROMAT μ SIC by DEKEMA

AUSTROMAT 624 i by DEKEMA


AUSTROMAT 674 by DEKEMA

AUSTROMAT 664 i by DEKEMA


AUSTROMAT 674 i by DEKEMA

AP ROTARY SYSTEMS

The AP Rotary Systems have been hand selected and tested by Master Dental Technician, Jörg Müller and a variety of other labs for specific needs and results.

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FOR A COMPREHENSIVE BUR LIST AND PRICES, PLEASE GO TO OUR WEBSHOP
WWW.APDENTAL.NET

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